



TEL:4006-871-227 Web:www.ybio.net Email:shybio@126.com

## 补体成分 1s(C1s)重组蛋白

适用生物 Homo sapiens (Human , 人)  
产品应用 SDS-PAGE; WB; ELISA; IP. 规格 50ug  
分子量 31.1kDa 价格 请咨询当地代理商  
纯度 > 95%  
来源 原核表达


### Recombinant Complement Component 1, S Subcomponent (C1s)

FOR IN VITRO AND RESEARCH USE ONLY, NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES!

Organism species	Homo sapiens (Human)
Product No.	YBB674Hu01
Source	Prokaryotic expression
Host	E. coli
Purity	> 95%
UOM	50ug
Predicted Molecular Mass	31.1kDa
Predicted isoelectric point	8.2
Applications	SDS-PAGE; WB; ELISA; IP.
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)



TEL:4006-871-227 Web:www.ybio.net Email:shybio@126.com



Subcellular Location	n/a
Residues	Ile438~Pro685 (Accession # P09871) with two N-terminal Tags, His-tag and T7-tag
Formulation	Supplied as lyophilized form in PBS, pH7.4, containing 5% trehalose, 0.01% sarcosyl

## SEQUENCES

The target protein is fused with two N-terminal Tags, His-tag and T7-tag, its sequence is listed below.

```
MGSSHHHHH SGLVPRGSH MASMTGGQM GRGSEF- IIG GSDADIKNFP WQVFFDNPWA G
GALINEYWV LTAHVVEGN REPTMYVGST SVQTSRLAKS KMLTPEHVFI HPGWKLEVP EGRTN
FDNDI ALVRLKDPVK MGPTVSPICL PGTSSDYNLM DGDGLISGW GRTEKRDRAV RLKAARLPV
A PLRKCKEVKV EKPTADAEAY VFTPNMICAG GEKGMDSCKG DSGGAFVAVQD PNDKTKFYAA G
LVSWGPPCG TYGLYTRVKN YVDWIMKTMQ ENSTP
```

## USAGE

Reconstitute in sterile PBS, pH7.2-pH7.4.

## STORAGE AND STABILITY

**Storage:** Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.



TEL:4006-871-227 Web:www.ybio.net Email:shybio@126.com

**Stability Test:** The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

### About the **MARKER (complimentary)**

**Effective Size Range:** 10kDa to 70kDa.

**Protein bands:** 10kDa, 14kDa, 18kDa, 22kDa, 26kDa, 33kDa, 44kDa and 70kDa.

**Double intensity bands:** The 26kDa, 18kDa, 10kDa bands are at double intensity to make location and size approximation of proteins of interest quick and easy.

**Ready-to-use:** No need to heat, dilute or add reducing agents before use.