



TEL:4006-871-227 Web:www.ybio.net Email:shybio@126.com

**YBB878Ca01 100μg**

**Recombinant Coagulation Factor VIII (F8)**

**Organism Species: Canis familiaris; Canine (Dog)**

***Instruction manual***

**FOR IN VITRO USE AND RESEARCH USE ONLY**

**NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES**

9th Edition (Revised in Jul, 2013)

## **[ PROPERTIES ]**

**Residues: Met2191~Cys2337 (Accession # O18806),  
with N-terminal His-Tag.**

**Host: *E. coli***

**Subcellular Location: Secreted, extracellular space.**

**Purity: >95%**

**Endotoxin Level: <1.0EU per 1μg (determined by the  
LAL method).**


**Formulation: Supplied as lyophilized form in PBS,  
pH7.4, containing 1mM DTT, 5% trehalose, 0.01%  
sarcosyl and preservative.**

**Predicted isoelectric point: 10.1**

**Predicted Molecular Mass: 18.2kDa**

**Applications: SDS-PAGE; WB; ELISA; IP.**

**(May be suitable for use in other assays to be determined by the end user.)**



## **[ USAGE ]**

**Reconstitute in sterile PBS, pH7.2-pH7.4.**



## **[ STORAGE AND STABILITY ]**

**Storage: Avoid repeated freeze/thaw cycles.**

**Store at 2-8°C for one month.**

**Aliquot and store at -80°C for 12 months.**

**Stability Test:** The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

## **[ SEQUENCES ]**

The target protein is fused with N-terminal His-Tag, its sequence is listed below.

MG HH HHH HSG SE F-ME SKA ISDA Q ITAS SYLSS M LATWSP SQ ARLH LQG RTNAW  
RPQANNPKIEW LQVDFRKTMK VTGITTQGVK SLLISMYVKE FLISSSQDGHNWTFLQNGK  
VKVFQGNRDS STPVRNRLEP PLVARYVRLH PQSWAHHIAL RLEVLGC

## **[ REFERENCES ]**

1. Cameron C., *et al.* (1998) *Thromb. Haemost.* 79:317-322.
2. Antonarakis SE. (1995) *Thromb. Haemost.* 74 (1): 322-8.
3. Lollar P. (2003) *Vox Sang.* 83. Suppl 1: 403-8.
4. Fay PJ. (2004) *Blood Rev.* 18 (1): 1-15.