

TEL:4006-871-227 Web:www.ybio.net Email:shybio@126.com

YBA133Po01 10µg

Recombinant Tumor Necrosis Factor Alpha (TNFa) Organism Species: Sus scrofa; Porcine (Pig) Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

kDa 94 66.2 45 33 26 20 14.4 15% SDS-PAGE

9th Edition (Revised in Jul, 2013)

[PROPERTIES]

Residues: Leu77^{Leu232} (Accession # P23563), with two N-terminal Tags, His-tag and T7-tag. Host: E. coli Subcellular Location: Cell membrane; Secreted Single-pass type II membrane protein. Purity: >95% Endotoxin Level: <1.0EU per 1µg (determined by the LAL method). Formulation: Supplied as lyophilized form in PBS, pH7.4, containing 5% sucrose, 0.01% sarcosyl. Predicted isoelectric point: 6.8 Predicted Molecular Mass: 21.0kDa Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)

[USAGE]

Bio11-227 TEL:4006-871-227

Web:www.ybio.net Email:shybio@126.com Reconstitute in sterile PBS, pH7.2-pH7.4.



[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37° C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCES]

The target protein is fused with two N-terminal Tags, His-tag and T7-tag, its sequence is listed below.

MG SSHH HHH H SSGLVPRG SH MAS MTG GQ QM GRG SEF- LR SS SQ TSDKP VAH VVANVKAEGQ LQWQSGYANA LLANGVKLKD NQLVVPTDGL YLIYSQVLFR GQGCPSTNVF LTHTISRIAV SYQTKVNLLS AIKSPCQRET PEGAEAKPWY EPIYLGGVFQ LEKDDRLSAE INLPDYLDFA ESGQVYFGII AL

[REFERENCES]

1. Drews R.T., et. al. (1990) J. Nucleic Acids Res. 18:5564-5564.

- 2. Kuhnert P., et. al. (1991) J.Gene 102:171-178.
- 3. Choi C.S., et.al. (1991) J.Anim. Biotechnol. 2:97-105.
- 4. Chardon P., et, al. (2001) J. Tissue Antigens 57:55-65.
- 5. Pauli U., et. al. (1989) J. Gene 81:185-191.