



TEL:4006-871-227 Web:www.ybio.net Email:shybio@126.com

YBC321Mu01 100 μ g

Recombinant Asporin (ASP_N)

Organism Species: *Mus musculus* (Mouse)

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

10th Edition (Revised in Jan, 2014)

[PROPERTIES]

Residues: Lys16~Lys373

Tags: Two N-terminal Tags, His-tag and GST-tag

Accession: Q99MQ4

Host: *E. coli*

Subcellular Location: Secreted, extracellular space,
extracellular matrix.

Purity: >95%

Endotoxin Level: <1.0EU per 1 μ g (determined by the
LAL method).


Formulation: Supplied as lyophilized form in PBS,
pH7.4, containing 5% trehalose, 0.01% sarcosyl.

Predicted isoelectric point: 8.8

Predicted Molecular Mass:

70.9kDa

Applications: SDS-PAGE; WB; ELISA; IP.





TEL:4006-871-227 Web:www.ybio.net Email:shybio@126.com

(May be suitable for use in other assays to be determined by the end user.)

[USAGE]

Reconstitute in sterile PBS, pH7.2-pH7.4.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCES]

The sequence of the target protein is listed below.

KPFFS PSHTALKNMM LKDMEDTDDD DNDDDDNSLF PTKEPVNPPF PFDLFPTCPF
GCQCYSRVVH CSDLGLTSVP NNIPFDTRMV DLQNNKIKEI KENDFKGLTS LYALILNNNK
LTKIHPKTF LTTKLRRLYL SHNQLSEIPL NLPKSLAELR IHDNKVKKI QKDTFKGMNAL
HVLEMSANPLENNGIEPGAF EGVTVFHIRI AEAKLTSIPK GLPPTLLELH LDFNKISTVE
LEDLKRYREL QRLGLGNRI TDIENGTAN IPRVREIHLE HNKLKKIPSG LQELKYLQII
FLHYNSIAKV GVNDFCPTVP KMKKSLYS AI SLFNNPMKYW EIQPATFRCV LGRMSVQLGN
VGK