



TEL:4006-871-227 Web:www.ybio.net Email:shybio@126.com

YBB329Mu01 100 $\mu$ g

Recombinant Microtubule Associated Protein 2 (MAP2)

Organism Species: *Mus musculus* (Mouse)

*Instruction manual*

FOR IN VITRO USE AND RESEARCH USE ONLY  
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

9th Edition (Revised in Jul, 2013)

## [ PROPERTIES ]

Residues: Thr1457~Ser1728 (Accession # P20357), with two N-terminal Tags, His-tag and GST-tag.

Host: *E. coli*

Subcellular Location: Cytoplasm, cytoskeleton.

Purity: >95%

Endotoxin Level: <1.0EU per 1 $\mu$ g (determined by the LAL method).


Formulation: Supplied as lyophilized form in PBS, pH7.4, containing 5% sucrose, 0.01% sarcosyl.

Predicted isoelectric point: 9.7

Predicted Molecular Mass: 61.2kDa

Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)





## [ USAGE ]

Reconstitute in sterile PBS, pH7.2-pH7.4.

## [ STORAGE AND STABILITY ]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

## [ SEQUENCES ]

The target protein is fused with two N-terminal Tags, His-tag and GST-tag, its sequence is listed below.

MSPILGYWKI KGLVQPTRLL LEYLEEKYEE HLYERDEGDK WRNKKFELGL EFPNLPYYID  
GDVKLTQSM A IIRYIADKHN MLGGCPKERA EISMLEGAVL DIRYGVSRIA YSKDFETLKV  
DFLSKLP EML KMFEDRLCHK TYLNGDHVTH PDFMLYDALD VVLYMDPMCL DAFPKLVCFK  
KRIEAI PQID KYLKSSKYIA WPLQGWA TFG GGDHPPKSD GSTSGSGHHH HHHSAGLVPR  
GSTAIGMKET AAKAFERQHM DSPDLGTLEV LFGQPLGSEF-TVSR DEVRRKKAVY  
KKAELAKKSE VQAHSPSRKL ILKPAIKYTR PTHLSCVKRK TTAASGDLAQ APGAFKQAKD  
KVTDGISKSP EKRSSLP RPS SILPPRRGVS GDREENSFSL NSSISSARRT TRSEPIRRAG  
KSGTSTPTTP GSTAITPGTP PSYSSRTPGT PGTPSYPRTP GTPKSGILVP SEKKVAIIRT  
PPKSPATPKQ LRLINQPLPD LKNVSKIGS TDNIKYQPKG GQVQIVTKKI DLSHVTSKCG  
SLKNIRHRPG GGRVKIES