

YBB389Hu01 50µg

# Recombinant Von Willebrand Factor A Domain Containing Protein 2 (vWA2)

#### Organism Species: Homo sapiens (Human)

#### Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

11th Edition (Revised in May, 2016)

# [PROPERTIES]

Source: Prokaryotic expression. Host: E. coli Residues: Ser457~Gly708 Tags: N-terminal His-Tag Tissue Specificity: Brain, Placenta. Subcellular Location: Secreted. **Purity: >98% Endotoxin Level:** <1.0EU per 1µg (determined by the LAL method). Traits: Freeze-dried powder Buffer formulation: 20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% sarcosyl, 5% Trehalose and Proclin300. Original Concentration: 200ug/mL Applications: SDS-PAGE; WB; ELISA; IP; CoIP; Reporter Assays; Purification; Amine Reactive Labeling. (May be suitable for use in other assays to be determined by the end user.) Predicted isoelectric point: 5.3 Predicted Molecular Mass: 30.5kDa Accurate Molecular Mass: 31kDa as determined by SDS-PAGE reducing conditions.



TEL:4006-871-227 Web:www.ybio.net Email:shybio@126.com

# [USAGE]

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

### [ STORAGE AND STABILITY ]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37<sub>o</sub>C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

## [SEQUENCE]

SEDE VAGPARHARA RELLLLGVGS EAVRAELEEI TGSPKHVMVY SDPODLFNOI PELOGKLCSR ORPGCRTOAL DLVFMLDTSA SVGPENFAOM OSFVRSCALO FEVNPDVTOV GLVVYGSOVO TAFGLDTKPT RAAMLRAISO APYLGGVGSA GTALLHIYDK VMTVORGARP GVPKAVVVLT GGRGAEDAAV PAOKLRNNGI SVLVVGVGPV LSEGLRRLAG PRDSLIHVAA YADLRYHODV LIEWLCGG

### [IDENTIFICATION]




Figure 1. SDS-PAGE

