

TEL:4006-871-227 Web:www.ybio.net Email:shybio@126.com

YB92127Hu0

1

Decorin (DCN)

Organism: Homo sapiens (Human)

Instruction manual

# FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN DIAGNOSTIC OR THERAPEUTIC PROCEDURES

4th Edition (Revised in February, 2012)

#### [ DESCRIPTION ]

**Protein Names:** Decorin **Synonyms:** DCN, SLRR1B

Species: Human

Size: 100µg

Source: Escherichia coli -derived Subcellular Location: Secreted.

### [PROPERTIES]

Residues: Asp31~Lys359 (Accession # P07585), with N-terminal His-Tag.

**Grade & Purity:** >97%, 37.8kDa as determined by SDS-PAGE reducing conditions. **Formulation:** Supplied as liquid form in Phosphate buffered saline(PBS), pH 7.4.

**Endotoxin Level:** <1.0 EU per 1µg (determined by the LAL method).

Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)

Predicted Molecular Mass: 37.8kDa Predicted isoelectric point: 8.61

### [PREPARATION]

Reconstitute in sterile PBS, pH7.2-pH7.4.



TEL: 4006-871-227 Web: www.ybio.net Email: shybio@126.com

### [STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

**Stability Test:** The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed.(Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

#### [SEQUENCES]

The target protein is fused with one N-terminal His-tag, its sequence is listed below. DEASGIGPEV PDDRDFEPSL GPVCPFRCQC HLRVVQCSDL GLDKVPKDLP PDTTLLDLQN NKITEIKDGD FKNLKNLHAL ILVNNKISKV SPGAFTPLVK LERLYLSKNQ LKELPEKMPK TLQELRAHEN EITKVRKVTF NGLNQMIVIE LGTNPLKSSG IENGAFQGMK KLSYIRIADT NITSIPQGLP PSLTELHLDG NKISRVDAAS LKGLNNLAKL GLSFNSISAV DNGSLANTPH LRELHLDNNK LTRVPGGLAE HKYIQVVYLH NNNISVVGSS DFCPPGHNTK KASYSGVSLF SNPVQYWEIQ PSTFRCVYVR SAIQLGNYK

## [REFERENCES]

- 1. Chen R., et al. (2009) J. Proteome Res. 8:651-661.
- 2. Bredrup C., et al. (2005) Invest. Ophthalmol. Vis. Sci. 46:420-426.
- 3. Danielson K.G., et al. (1993) Genomics 15:146-160.
- 4. Fisher L.W., et al. (1987) J. Biol. Chem. 262:9702-9708.
- 5. Krusius T., et al. (1986) Proc. Natl. Acad. Sci. U.S.A. 83:7683-7687.