



TEL:4006-871-227 Web:www.ybio.net Email:shybio@126.com

YBC756Hu01 100µg

Recombinant Prosaposin (PSAP)

Organism Species: Homo sapiens (Human)

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

10th Edition (Revised in Jan, 2014)

[PROPERTIES]

Residues: Gly17~Asn524

Tags: Two N-terminal Tags, His-tag and T7-tag

Accession: P07602

Host: *E. coli*

Subcellular Location: Lysosome.

Purity: >95%

Endotoxin Level: <1.0EU per 1µg (determined by the LAL method).


Formulation: Supplied as lyophilized form in PBS, pH7.4, containing 5% trehalose, 0.01% sarcosyl.

Predicted isoelectric point: 5.1

Predicted Molecular Mass: 60.2kDa

Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)



[USAGE]

Reconstitute in sterile PBS, pH7.2-pH7.4.



[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCES]

The sequence of the target protein is listed below.

GPVL GLKECTRGS A VWCQNVKTAS DCGAVKHCLQ TVWNKPTVKS LPCDICKDVV
TAAGDMLKDN ATEEEILVYL EKTCDWLPKP NMSASCKEIV DSYLPVILDI IKGEMSRPGE
VCSALNLCES LQKHLAELNH QKQLESNKIP ELDMTEVVAP FMANIPLLLY PQDGPRSKPQ
PKDNGDVCQD CIQMVTDIQT AVRTNSTFVQ ALVEHVKEEC DRLGPGMADI CKNYISQYSE
IAIQMMMHHMQ PKEICALVGF CDEVKEMPMQ TLVPAKVASK NVIPALELVE PIKKHEVPAK
SDVYCEVCEF LVKEVTKLID NNKTEKEILD AFDKMCSKLP KSLSEECQEV VDTYGSSILS
ILLEEVPEL VCSMLHLCSG TRLPALTVHV TQPKDGGFCE VCKKLVGYLD RNLEKNSTKQ
EILAALEKGC SFLPDYQKQ CDQFVAEYEP VLIEILVEVM DPSFVCLKIG ACPSAHKPLL
GTEKCIWGPS YWCQNTETAA QCNAVEHCKR HVWN