TEL:4006-871-227 Web:www.ybio.net Email:shybio@126.com

YB93950Mu01 100µg

Pregnane X Receptor (PXR)

Organism: Mus musculus (Mouse)

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

7th Edition (Revised in May, 2013)

## [ PROPERTIES ]

Residues: Tyr205-Gly390 (Accession # 054915), with two Nterminal Tags, His-tag and GST-tag.

Host: E. coli

Subcellular Location: Nucleus.

Purity: >95%

Endotoxin Level: <1.0EU per 1μg (determined by the LAL

method).


Formulation: Supplied as lyophilized form in PBS, pH7.4,

containing 5% sucrose, 0.01% sarcosyl.

Predicted isoelectric point: 6.6 Predicted Molecular Mass: 48.7kDa

Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)





TEL:4006-871-227 Web:www.ybio.net Email:shybio@126.com

## [ USAGE ]

Reconstitute in sterile PBS, pH7. 2-pH7. 4.

## [ STORAGE AND STABILITY ]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at  $-80^{\circ}$ C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at  $37^{\circ}\mathrm{C}$  for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

## [ SEQUENCES ]

The target protein is fused with two N-terminal Tags, His-tag and GSTtag, its sequence is listed below.

MRNKKFELGL EFPNLPYYID GDVKLTQSMA IIRYIADKHN MLGGCPKERA EISMLEGAVI.

DIRYGVSRIA YSKDFETLKV DFLSKLPEML KMFEDRLCHK TYLNGDHVTH PDFMLYDALD

VVLYMDPMCL DAFPKLVCFK KRIEAIPQID KYLKSSKYIA WPLQGWQATF GGGDHPPKSD

G S T S G S G H H H H H H H S A G LV P R G S TA I G M K E T A A A K F E R Q H M D S P D L G T L E V

LFQGPLGSEF-YMFKGV INFAKVISYF RDLPIEDQIS LLKGATFEMC ILRFNTMFDT

ETGTWECGRL AYCFEDPNGG FQKLLLDPLM KFHCMLKKLQ LHKEEYVLMQ AISLFSPDRP

GVVQRSVVDQ LQERFALTLK AYIECSRPYP AHRFLFLKIM AVLTELRSIN AQQTQQLLRI

QDSHPFATPL MQELFSSTDG