



TEL:4006-871-227 Web:www.ybio.net Email:shybio@126.com

YB90004Mu01

Angiotensin I Converting Enzyme (ACE)

Organism: *Mus musculus* (Mouse)

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY

NOT FOR USE IN DIAGNOSTIC OR THERAPEUTIC PROCEDURES

4th Edition (Revised in August, 2012)

[DESCRIPTION]

Protein Names: Angiotensin I Converting Enzyme

Synonyms: ACE, Dcp1

Species: Mouse

Size: 100 μ g

Source: *Escherichia coli*-derived

Subcellular Location: Secreted, Cell membrane; Single-pass type I membrane protein.

[PROPERTIES]

Residues: Tyr504~Asp597 (Accession # P09470), with two N-terminal Tags, His-tag and GST-tag.

Grade & Purity: >95%, 38 kDa as determined by SDS-PAGE reducing conditions.

Formulation: Supplied as liquid form in Phosphate buffered saline(PBS), pH 7.4.

Endotoxin Level: <1.0 EU per 1 μ g (determined by the LAL method).

Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)

Predicted Molecular Mass: 37.6 kDa

Predicted isoelectric point: 6.8

[PREPARATION]

Reconstitute in sterile PBS, pH7.2~pH7.4.



TEL:4006-871-227 Web:www.ybio.net Email:shybio@126.com

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCES]

The target protein is fused with two N-terminal Tags, His-tag and GST-tag, its sequence is listed below.

MRNKKFELGLEFPNLPYYIDGDVKLTQSMAIIRYIADKHNMLGGCPKERAEISMLEGAVLDIRYGVSR IAYSKDFETLKVDFLSKLP EML

KMFEDRLCHKTYLNGDHVTHPDFMLYDALDVVLYMDPMCLDAFPKLVCFKKRIEAI PQIDKYLKSSKYIAWPLQGWQATFGGGDHPP

KSDGSTSGSGHHHHHSAGLVPRGSTAIGMKETA AAKFERQHMDSPDLGTLEVL FQGPLGSEF-YQGICPP

VARNETHFDA

GAKFHIPNVT PYIRYFVSFV LQFQHQALC KEAGHQPLH QCDIYQSAQA GAKLKQVLQA GCSR PWQEV L KDLVGSD